



Mercury Series

Portable Oxygen Concentrator

OX-2A



Dual mode Continuous / Pulse Dose

Battery Duration

Oxygen Concentration	87%-96% at all settings
Flow Settings and Pulse Volumes	Pulse Mode 1 -12 settings (Average of 20 consecutive breaths)
	Continuous Mode 0.5 -2 Liters per mins (3 minute running average) Note: Max recommended flow: 2 LPM
User Interface	Power switch, Push buttons, back-light liquid crystal display (LCD)
Dimensions	200 (L) X 230 (W) X 325 (H) mm
Weight	5.4 kg / 4 kg (Elite)
Sound Level	46 dBA max. at continuous Mode Setting of 2 liter Front 1M
Audio Alarm	50 decibels (nominal) at 39 in. (1 m)
Outlet Pressure	8.5 psi ± 0.5 psi

Setting	Pulse Dose (20 BPM)		Continuous	
	ML	Mins	L/Min	Mins
1	6	300	0.5	240
2	12	235	1	165
3	18	235	1.5	120
4	24	200	2	75
5	30	200		
6	36	165		
7	42	165		
8	48	145		
9	54	135		
10	60	120		
11	66	115		
12	72	110		

*Based on an atmospheric pressure of 14.7 psi (101 kPa) at 70°F (21°C)

