

INSPIRATION™

Revolutionary Design:

The Labex Inspiration Electrolarynx is the most advanced electrolarynx device on the market.

It combines the latest technology with intelligent design.

Our groundbreaking Touch Technology allows you to change the pitch of your voice during speaking – simply slide a finger on the control panel!

This permits the user to express emotions by changing intonation.

A detailed light indication makes the device easy to adapt to, while our Smart Start feature lets the electrolarynx to start from the moment it touches the neck.

- The powerful built-in rechargeable battery allows for 24 hours of talk time
- Full charge battery time is 2.5 hours with a factory-kit charger
- The solid polished stainless body is easy to use and clean
- Superior voice-matching capabilities
- Easy-to-read LED indicator
- Settings are intuitive — easy to adapt to and use!



HARMONY™

Thanks to our simple two-button design, each button can be configured to different volume level and tone.

It is useful to have a button for outdoor use in a noisy environment. Second button could be set to a lower volume to allow you communicate indoor or in situations when you want to speak quietly.

- Realistic and clear speech
- High quality Industrial Design
- Strong metal case
- One charge for more than 20 hours speech time
- Well felt buttons
- Built-in lithium-ion battery high throughput and long service life.



Digital

Comfort

Digital Electolaynxes

Labex Digital and Labex Comfort –

Voice-forming devices

Labex Digital and Labex Comfort -

some of the most compact and light voice-forming devices on the market.

They weigh only 98 grams, thereby providing maximum mobility!

Labex Comfort and Labex Digital work standard 9V batteries.

Included

two 9V Li-Ion rechargeable rechargeable batteries and charger device.

Comfort and Digital models made from durable, shockproof plastic.

Thanks to our simple design, each button can easily be tuned to your level volume and tone..



- Two speech buttons for Comfort.
- One speech button for Digital.
- Each button has its own values of volume and pitch.
- Digital control allows talk for 4 hours per battery charge.
- Convenient two-stage warning when the battery is low.
- IP54 degree of moisture protection and dustproof
- Very rugged case.
- Dimensions 30 x 30 x 103 mm
- Weight 72 g without, and 98 g with battery