

Oxygen concentrator

OX 5A/ OX 10A

OXTm



- ✓ High flow rate-1-10 Liters
- ✓ Purity for all flow at 93±3%
- ✓ The most reliable, efficient, and convenient source of supplemental oxygen
- ✓ Provide supplemental low flow oxygen therapy in the home and patient care facilities

	OX-5A	OX-10A
Delivery Rate (LPM)	0.5-5	2~10
Maximum Flow (LPM)	5	10
Output Pressure (KPa)	58 ± 6 (8.5± 0.5 Psi)	80 ± 8 (12 ± 1 Psi)
Electrical Requirement (Voltage/Hz)	110VAC , 220VAC 50/60Hz	110VAC , 220VAC 50/60Hz
Purity (%)	93±3(0.5-5.0 L)	93±3 (2-9 L) / 90±3 (10 L)
Power Consumption (Watt)	300 (Aver.)	580 (Aver.)
Weight (Kg)	15.2	19.6
Noise Level DB(A)	Front 1 M < 42	Front 1 M < 50
Dimensions (Mm)	368 (L) X 342(W) X 572 (H)	368 (L) X 342 (W) X 572 (H)
OCI (Oxygen Concentration Indicator)	YELLOW: 82% Below RED: 70% Below	YELLOW: 82% Below RED: 70% Below
Flow meter	Float	Float
Nozzle	1	1
Warranty	3 Years	1 Years